
Pinnacle Series

Pinnacle-Serie

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English – page 3

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Contents

Warning and Safety Information . . .	4	Troubleshooting Guide	18
Specifications	5	Balance Operating Menu	19
Getting Started	8	Accessories	22
Calibration and Adjustment	12	Declaration of Conformity	23
Operation	14	Confirmation of an EC Type-Approval Certificate	25
Counting	17	Plates and Markings	26
Target weighing/ Weighing in Percent	17		
Checkweighing	17		
Animal Weighing	17		
Formulation	17		
Statistics	17		
GLP/ISO	17		
Printing	17		

Denver Instrument Pinnacle Series balances meet high quality standards, and with proper care and treatment will provide many years of trouble-free service.

Please read this manual carefully before operating your new Denver Instrument balance.

Warning and Safety Information

Calibrate (adjust) your balance using a reference weight of the appropriate accuracy class. A measuring instrument is only as accurate as the standard that was used as a reference in calibration. For assistance in the selection of reference weights, please contact Denver Instrument Company.

Caution

The operator shall be responsible for any modifications to Denver Instrument equipment and for any connections of cables or equipment not supplied by the Denver Instrument Company and must check and, if necessary, correct these modifications and connections.

Caution

Use of this product in a manner other than that specified by the manufacturer may impair any safety protection provided by the equipment.

Warning

Do not drop objects (such as weights or samples) on the weighing pan, as this will damage the balance's internal mechanisms.

Warning

Never lift balance by the weighing pan or pan support as this may cause damage to internal mechanisms.

Always lift and transport the balance by its base, including when removing the balance from the packaging.

Note:

Loads that exceed the balance's maximum capacity by 15% can damage the weighing system.

Note:

This product may contain a rechargeable battery.

The user is responsible for disposing of the equipment, including the battery, in accordance with local ordinances and regulations.

Specifications

Standard Models

Precision and Analytical Balances

Models	P-114*	P-214*	P-314*
Capacity	110 g	210 g	310 g
Readability	0.1 mg	0.1 mg	0.1 mg
Taring range	0 to 110 g	0 to 210 g	0 to 310 g
Repeatability, (s)	0.1 mg	0.1 mg	0.1 mg
Linearity	0.2 mg	0.2 mg	0.3 mg
Response time	3 seconds	3 seconds	3 seconds
Pan size	79mm Ø (3.1 in.)	79mm Ø (3.1 in.)	79mm Ø (3.1 in.)

Precision Balances

Models	P-203*	P-403*	P-603D*	P-402*	P-602*	P-2002*
Capacity	200 g	400 g	100/600 g	400 g	600 g	2000 g
Readability	0.001 g	0.001 g	0.001/0.01g	0.01g	0.01g	0.01 g
Taring range	0 to 200 g	0 to 400 g	0 to 100/600 g	0 to 400 g	0 to 600 g	0 to 2000g
Repeatability, (s)	0.001 g	0.001 g	0.002/0.01 g	0.01 g	0.01 g	0.01 g
Linearity	0.002 g	0.002 g	0.003/0.02 g	0.02 g	0.02 g	0.02 g
Response time	3 seconds	3 seconds	4 seconds	3 seconds	3 seconds	3 seconds
Pan size	114mm (4.5 in.)	Ø114mm Ø (4.5 in.)	114mm Ø (4.5 in.)	114mm Ø (4.5 in.)	114mm Ø (4.5 in.)	178 x 178mm (7.0 x 7.0 in.)

Models	P-4002*	P-4002D*	P-8002D	P-6001*	P-8001
Capacity	4000 g	400/4000 g	800/8000 g	6000 g	8000 g
Readability	0.01 g	0.01/0.1 g	0.01/0.1 g	0.1 g	0.1 g
Taring range	0 to 4000 g	0 to 400/4000 g	0 to 800/8000 g	0 to 6000 g	0 to 8000 g
Repeatability, (s)	0.01 g	0.02/0.1 g	0.02/0.1 g	0.1 g	0.1 g
Linearity	0.02 g	0.03/0.2 g	0.03/0.2 g	0.2 g	0.2 g
Response time	3 seconds	4 seconds	4 seconds	3 seconds	3 seconds
Pan size	178 x 178mm (7.0 x 7.0 in.)	178 x 178mm (7.0 x 7.0 in.)	178 x 178mm (7.0 x 7.0 in.)	178 x 178mm (7.0 x 7.0 in.)	178 x 178mm (7.0 x 7.0 in.)

* PI models with internal calibration

Specifications

Verified Models with an EC Type-Approval Certificate

Precision Analytical Balances

Models	PI-114A	PI-214A	PI-314A	PI-225DA
Accuracy class ¹⁾	I	I	I	I
Type	BD BP 100	BD BP 100	BD BP 100	BC BP 100
Max. capacity, Max ¹⁾	110 g	210 g	310 g	60/220 g
Scale interval, d ¹⁾	0.1 mg	0.1 mg	0.1 mg	0.01/0.1 mg
Taring range (subtractive)	0 to 110 g	0 to 210 g	0 to 310 g	0 to 220 g
Verification scale interval, e ¹⁾	1 mg	1 mg	1 mg	1 mg
Minimum capacity, Min ¹⁾	10 mg	10 mg	10 mg	1 mg
Response time	3 seconds	3 seconds	3 seconds	6/3 seconds
Pan size	79 mm Ø (3.1 in.)	79 mm Ø (3.1 in.)	79 mm Ø (3.1 in.)	79 mm Ø (3.1 in.)

Precision Balances

Models	PI-203A	PI-403A	PI-602A	PI-2002A	PI-4002A	PI-6001A
Accuracy class ¹⁾	II	II	II	II	II	II
Type	BD BP 200	BD BP 200	BD BP 200	BD BP 200	BD BP 200	BD BP 200
Max. capacity, Max ¹⁾	200 g	400 g	600 g	2000 g	4000 g	6000 g
Scale interval, d ¹⁾	0.001 g	0.001 g	0.01 g	0.01 g	0.01 g	0.1 g
Taring range (subtractive)	0 to 200 g	0 to 400 g	0 to 600 g	0 to 2000g	0 to 4000 g	0 to 6000 g
Verification scale interval, e ¹⁾	0.01 g	0.01 g	0.1 g	0.1 g	0.1 g	1 g
Minimum capacity, Min ¹⁾	0.02 g	0.02 g	0.5 g	0.5 g	0.5 g	5 g
Response time	3 seconds	3 seconds	3 seconds	3 seconds	3 seconds	3 seconds
Pan size	114mm Ø (4.5 in.)	114mm Ø (4.5 in.)	114mm Ø (4.5 in.)	178x178mm (7.0 x 7.0 in.)	178x178mm (7.0 x 7.0 in.)	178x178mm (7.0 x 7.0 in.)

¹⁾ Council Directive 90/384/EEC on non-automatic weighing instruments used within the European Economic Area

General Specifications:

Dimensions (WxDxH), analytical models:	371 x 239 x 338mm (14.6 x 9.4 x 13.3 in.)
Dimensions (WxDxH), precision models:	371 x 239 x 97mm (14.6 x 9.4 x 3.8 in.)
Weighing chamber dimensions:	211 x 203 x 249mm (8.3 x 8.0 x 9.8 in.)
Allowable operating temperature range:	10° - 30°C (50° - 86°F)
Storage temperature:	-10° - 30°C (14° - 86°F)
Humidity:	< 90% RH
Net weight (analytical balances):	6.8 kg (15 lbs)
Net weight (precision balances):	4.54 kg (10 lbs)
Power requirements/frequency: AC:	120-240V 50/60Hz; other voltages available

NOTE:

The battery is not a user-serviceable part. If the battery does not recharge while the AC adapter is plugged in, contact an authorized service technician.

Caution:

Use only the AC adapter/charger supplied with balance.

Contact Denver Instrument for replacement.

Getting Started

Installation Instructions

Your balance is designed to provide reliable weighing results under normal ambient conditions. When choosing a location to set up your new balance, observe the following so that you will be able to work with added speed and accuracy:

- Set up the balance on a even, stable surface.
- Avoid exposing the balance to extreme heat radiation from heaters or direct sunlight
- Due not use the balance at temperatures exceeding the allowable operating temperature range (+10° to +30°C; 50° to 86°F) as this could negatively affect balance operation and accuracy.
- Protect the balance from direct exposure to drafts.
- Protect the balance from aggressive chemical vapors.
- Avoid strong magnetic fields generated by other devices.
- Avoid exposing the balance to extreme vibrations during weighing
- Line voltage to the balance should be reasonably constant and free from fluctuation.

Conditioning the Balance

Moisture in the air can condense on the surfaces of a cold balance whenever it is brought to a substantially warmer place. If you transfer the balance to a warmer area, make sure to condition it for about 2 hours at the new ambient temperature, leaving it unplugged from the power supply. Afterwards, keep the balance connected to the power supply.

Warranty


In the unlikely event that your balance should require servicing within the 3-year warranty period, please contact the office listed below for your country, and state the model number, serial number and the problem:
Federal Republic of Germany:

DENVER INSTRUMENT GmbH
Robert-Bosch-Breite 10
37079 Goettingen
Germany
Telephone: +49 (0) 551/209 773-0
Telefax: +49 (0) 551/209 773-9
Internet:
<http://www.DenverInstrument.com>
E-mail: info@DenverInstrument.de

Seal on Balances Verified for Use in Legal Metrology in the EU*:

A protective seal is required for verified balances (see page 24). The protective seal consists of a sticker. This seal will be irreparably damaged if you attempt to remove it. If the seal is broken, the validity of verification will become void and you must have your balance re-verified.


* including the Signatories of the Agreement on the European Economic Area

Analytical Balances and Draft Shield

To attach draft shield, place it upon the balance so the grooves line up when the draft shield is approximately at a 30° angle. Pull the pin at the rear of the draft shield and turn clockwise until it is in line with the balance. Release the pin. Place the components listed below inside the chamber in the order given:

- Stainless steel platform
- Spill ring
- Weighing pan



Connecting Model PI-225DA to the Electronics Box

- Plug the male connector on the cable into the female connector on the electronics box

Precision Balances with Round Pan

Place the components listed below on the balance in the order given:

- Shield ring
- Base pan
- Weighing pan
- Draft shield (optional)
- Draft shield cover (optional)

Place the weighing pan on the balance so the rubber feet line up with the holes

Precision Balances with Square Pan


- Weighing pan

Connecting the Balance to AC Power

After conditioning the balance (see page 8), simply plug the AC adapter into the rear of the balance and connect the adapter to an appropriate power outlet; the balance switches on automatically. To avoid extended warm-up periods, leave the balance plugged in and switched on at all times.

PI-225DA

- Insert the right-angle plug from the AC adapter into the jack on the electronics box.

Back panel

Leveling the Balance

Level the balance any time it is installed in a new location. Make sure all four feet touch the work surface. Adjust the leveling feet until the air bubble is centered in the level indicator.

The number of leveling feet depends on the balance model:

Analytical models: 2 front feet

Precision models with round pan: 2 front feet

Precision models with square pan: 4 feet, one in each corner

* = not on verified PI ... A models

Calibration and Adjustment

Calibration technically means to determine the difference between the weighing instrument readout and the actual weight on the platform to determine the accuracy. Adjustment means to bring a weighing instrument to the level of accuracy required for its use and adjusted.

Although the term “calibration” is used throughout this manual, as is common in the USA, we are referring to “adjustment” (see first paragraph).

Your balance was calibrated at the factory; however, it must be re-calibrated upon installation and subsequently on a regular basis. The recommended calibration interval is once per week, using a certified calibration weight. Calibration must also be performed any time the balance has been moved, leveled, or disconnected from power, as well as following a significant change in ambient temperature. Additional calibration may be required for compliance with ISO/GLP standards.

Note on Verified Balances Used in Legal Metrology:

Before using your balance in legal metrology, perform internal calibration at the place of installation.

External Calibration/Adjustment*


1. Remove all items from the balance.
2. Press **TARE**.
3. Carefully center the calibration weight on the weighing pan.
4. Press the **⏻** key.
5. The display shows “**CAL**” and the weight.

* on verified PI ... A models, external calibration must be performed by a Denver service technician.

6. When calibration is complete, the weight value is displayed and the balance returns to the normal weighing mode.

Internal Calibration (PI models only)

1. Remove all items from the balance.
2. Press **TARE**.
3. Press the **⏻** key.
4. The display shows “CAL Internal” and calibration is performed automatically.

5. When calibration/adjustment is complete, the balance automatically returns to the normal weighing mode.


Note:

If “Internal Cal” is deactivated (“Off” in the operating menu), only external calibration can be performed.

If calibration/adjustment cannot be completed, an error message is displayed. See “Troubleshooting Guide” for details on error codes and messages.


DenverCal: Automatic Calibration/ Adjustment of PI...A Models

"DenverCal" flashes on the display to indicate that the ambient temperature has changed, or the specified time interval has elapsed, since the last calibration. The balance is ready to perform automatic calibration/adjustment.


Automatic calibration/adjustment is performed once the following two criteria have been met:

The load has not changed within the last 2 minutes. The "DenverCal" signal has been flashing for at least 15 seconds

"DenverCal" continues to flash until internal calibration/adjustment has been completed successfully.

Operation

Warm-up Time

To deliver exact results, the balance must be connected to AC power for at least 2 hours before weighing. Only after this time will the balance have reached the required operating temperature. We recommend leaving balance plugged in so that all components are at operating temperature and the balance is ready for use at any time.

Note on Verified Balances Used in Legal Metrology (PI...A Models):

- The balance must warm up for at least 24 hours after initial connection to AC power.
- The type-approval certificate for verification applies to non-automatic weighing instruments only; for automatic operation with or without auxiliary measuring devices, you must comply with the regulations applicable in your country.
- The temperature range indicated on the verification ID label must not be exceeded during operation.


Example:

BD BP 200

Ⓜ 0 to +40°C

Taring the Balance

1. Place a container on the weighing pan. The balance registers the weight of the container.
2. Press the TARE key.
3. "TARING" is displayed to indicate that the balance is being tared.
4. The display shows 0.0000 grams (or selected weight unit and resolution) once taring is concluded successfully.
5. When the readout is stable, the Denver Instrument logo appears in the upper left corner of the display. When unstable, a "U" is displayed.

The following features/applications are available in your balance:

- Calibration/Adjustment
- Target/Weighing/Weighing in Percent
- Animal Weighing
- GLP/ISO
- Counting
- Checkweighing
- Statistics
- Formulation

Accessing the operating menu:

1. Press the Menu key .
2. Press the ▼ key until the desired function or application is highlighted.
3. Press ► or the  key.
4. Press ▼ until the desired function is highlighted.
5. Use the arrow keys to select the desired menu item(s).
6. Press the  key to confirm input.

Note:

Press the TARE key to exit the menu at any point without saving changes.

Alphanumeric input:

1. Use the arrow keys to select the desired input field.
2. With the cursor in the input field, press ▼ to go to the first letter in alphabet, or ▲ to scroll through numbers.
3. Press the ► key to go to the next character(s).
4. Repeat this process until the desired value is entered.
5. Press the  key to confirm input and return to the main menu.

Press the ◀ key to return to the previous character. You can press and hold arrow keys to scroll quickly.

Counting

With the Counting application, you can determine a number of parts of approximately equal weight. Place the reference sample quantity of parts on the balance to determine the reference weight first; the next load on the balance is then calculated as a number of pieces having the same weight, based on this reference value. The readout shows both the piece count and the corresponding weight, in the primary weight unit.

Target Weighing/Weighing in Percent

This application allows you to obtain weight readouts in a percentage proportional to a reference weight. The balance will display both the percentage and the total weight in primary weight unit.

Checkweighing

This application gives you plain-language responses indicating the relation of the load on the balance to user-defined target and tolerance limits.

Animal Weighing

The Animal Weighing application is used for measuring animals that move continuously during weighing. This feature can also be used for measurements taken in environments with extreme vibrations and/or drafts.

"L" on the display indicates that the result is shown as a locked-in display value.

Formulation

This feature allows you to weigh out several different components in a single container. Simply place your container on the pan and press the TARE key. Add components one at a time, pressing ENTER after each component. The weight of each individual component is displayed as the primary weight, and the total weight is shown in the secondary weight display.

Statistics

Your balance is equipped with internal memory for 250 data records. You can calculate and print various statistics, including number of values, minimum weight, maximum weight, range, average, standard deviation and total weight.

GLP/ISO

Your balance has a number of features that let you configure data records for compliance with GLP and ISO documentation requirements. Activate the GLP/ISO header to include all required data on printouts.

Printing

The balance has a bi-directional RS-232 port for communication with serial peripheral devices, such as a printer or computer. In addition, the balance excluding PI... A models has a USB port* not on verified PI... A models for connection to a computer.

Ambient Conditions

You can adapt the balance by selecting the appropriate menu setting, to compensate for ambient conditions such as vibration and drafts.

* = not on verified PI... A models

System

The System menu lets you configure various features of your balance. You can assign password protection for menu settings, change the acoustic signal, restore factory defaults, and adjust the RS-232 interface settings for use with a computer or printer.

Restoring Factory Defaults


To restore all settings to the factory defaults:

1. Press the **Menu** key.
2. To open the System menu, press ▼ repeatedly until “System” is displayed, then press ►.
3. To restore the factory settings in the menu, press ▼ repeatedly until “Defaults” is displayed, then press ►.

Troubleshooting Guide


OVER (CAPACITY ERROR)	The load exceeds the balance capacity.	Unload the balance or remove obstruction.
	Display capacity is exceeded.	Decrease weight on balance.
UNDER	The load is too light.	Check pan position. Unload the balance and switch it off and then on again.
CAL ADD WEIGHT	Calibration/adjustment started without weight on pan.	Add weight and press CAL key again.
CAL OUT OF RANGE	Weight cannot be detected.	Unload pan and add appropriate weight.
Other errors	Error has occurred.	Switch the balance off and then on again. If the error persists, contact a service technician.

Balance Operating Menu


* = not on verified PI...A models

Balance Operating Menu

* = not on verified PL... A models

¹) = PL... A models: only when readout is stable

* = setting cannot be changed on PI...A models

Accessories

Analytical draft shield chamber for milligram models	602643.1
In-use cover (analytical models and models with round pan)	602619.1
In-use cover (square pan models)	602620.1
Printer	901042.1
Cable with 9-pin connector for balance	902225.1
Weigh below hanger*	77000440.8

Calibration weights: contact your distributor for a complete list

Complete operating instructions available at:
www.denverinstrument.com

* = not allowed in applications involving legal metrology
(such as legal for trade)

CE Declaration of Conformity to Council Directives 89/336/EEC and 73/23/EEC

**The electronic precision weighing
instrument of the series**

P/PI-.....

meets the requirements of the test standards listed below, in conjunction with the associated power supplies, auxiliary peripheral devices and installation equipment listed in Annex A2 (see Annex A1 for a technical description and variants).

1. Electromagnetic Compatibility

1.1 Source for 89/336/EEC: EC Official Journal, No. 2001/C105/03

EN 61326-1 Electrical equipment for measurement, control and laboratory use
EMC requirements
Part 1: General requirements

Limitation of emissions: Industrial areas, Class

A

**Defined immunity to interference: Minimum
immunity test requirements, continuous
unmonitored operation.**

2. Safety of Electrical Equipment

2.1 Source for 73/23/EEC: EC Official Journal, No. 2001/C106/03

EN 61010 Safety requirements for electrical equipment for measurement, control and laboratory use
Part 1: General requirements

Denver Instrument GmbH
Robert-Bosch-Strasse 10
37079 Göttingen, Germany
2003

Declaration of Type Conformity to Directive No. 90/384/EEC

The non-automatic weighing instruments with the EC type-approval certificate no.:

D03-09-016

Model	Type	Accuracy Class
PI-225DA	BC BP 100	(I)
PI-314A	BD BP 100	(I)
PI-214A	BD BP 100	(I)
PI-114A	BD BP 100	(I)
PI-403A	BD BP 200	(II)
PI-203A	BD BP 200	(II)
PI-602A	BD BP 200	(II)
PI-4002A	BD BP 200	(II)
PI-2002A	BD BP 200	(II)
PI-6001A	BD BP 200	(II)

Correspond to the production model described in the EC type-approval certificate and to the requirements of the Council Directive 90/384/EEC as amended.

The weighing instruments were initially verified under number 0111 by the manufacturer at the place of manufacture.

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Confirmation of an EC type-approval certificate

für die Firma, / for the company
DENVER INSTRUMENT GmbH, Robert-Bosch-Breite 10, 37079 Goettingen, Germany.

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

Siegel
Seal


Plates and Markings

Example of a descriptive plate on a weighing instrument already verified

Example of a plate with model designation

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

PI110603d

Typ: BD BP 100, BD BP 200
EG Bauartzulassung D03-09-016

Copyright by Denver Instrument GmbH, Goettingen, Germany.
All rights reserved. No part of this publication may be reprinted
or translated in any form or by any means without the prior
written permission of Denver Instrument GmbH.

The status of the information, specifications and illustrations in
this manual is indicated by the date given below. Denver Instrument
GmbH reserves the right to make changes to the technology, fea-
tures, specifications, and design of the equipment without notice.

Status: August 2004, Denver Instrument GmbH, Goettingen,
Germany.

Europe, Asia and Australia:

**Denver Instrument GmbH
Robert-Bosch-Breite 10
37079 Goettingen, Germany
Tel: +49-551-20977-30
Fax: +49-551-20977-39**

Internet:

www.denverinstrument.com

U.K. and Ireland:

**Denver Instrument Company
Denver House, Sovereign Way
Trafalgar Business Park
Downham Market
Norfolk PE38 9SW England
Tel: +44-136-63862-42
Fax: +44-136-63862-04**

North and South America:

**Denver Instrument Company
6542 Fig Street
Arvada, Colorado 80004
1-800-321-1135
Tel: +1-303-431-7255
Fax: +1-303-423-4831**